

Anti SARS-CoV-2(2019-nCoV) RBD Neutralizing antibody(LY-CoV555)**Basic information:**

Catalog No.:	UDA9009-5	Source:	CHO cells
Clonality:	Monoclonal	Isotype:	Human IgG

Useful Information:

Applications: WB, ELISA

Specificity: Recognizes SARS-CoV-2 RBD Protein

Form: Liquid

Buffer: PBS, pH7.5

Immunogen: Recombinant SARS-CoV-2 RBD protein

Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles.
Store at 2 to 8°C for one week .
Store at -20 to -80°C for twelve months from the date of receipt.

Note: Test systems and platform change may cause different results, please contact us if you have any questions.

This product is sold for research use only. Standard Laboratory Practices should be followed when handling this material.