

Recombinant SARS-CoV-2 (2019-nCoV) S Protein RBD(K417N, E484K, N501Y), His Tag


Basic information:

Catalog No.:	UDP9009-2	Expression system:	Mammalian cells
Form:	Liquid	Buffer:	PBS, pH7.5
Purity:	>90% as determined by SDS-PAGE		
Actual Molecular Mass:	35kDa		

Useful Information:

Applications: Immunogen
Description: Recombinant SARS-CoV-2 S Protein RBD(K417N, E484K, N501Y) is produced by Mammalian cells expression system and the target gene is expressed with C-His Tag

SDS-PAGE image:



Storage: Use a manual defrost freezer and avoid repeated freeze thaw cycles.
 Store at 2 to 8°C for one week .
 Store at -20 to -80°C for twelve months from the date of receipt.

Note: Test systems and platform change may cause different results, please contact us if you have any questions.
 This product is sold for research use only. Standard Laboratory Practices should be followed when handling this material.